

Anti-MgII (MagI) Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q99685

Item Number: CM0006942



Introduction

High-specificity rabbit monoclonal antibody targeting human MGLL (MAGL, HU-K5). Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism and endocannabinoid signaling research.

[Learn More](#)

Parameter	Value
Product Name	MGLL Rabbit mAb
Remarks/Alias	MGL, HUK5, MAGL, HU-K5, MGLL
Species	Human
GeneID (Human)	11343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MGLL
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99685
Protein Weight	33kDa
Shipping	Ice bag