

Anti-Mgll Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99685

Item Number: CM0029800



Introduction

Rabbit polyclonal antibody targeting MGLL, suitable for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Explore monoglyceride lipase function in endocannabinoid signaling, lipid metabolism, and cancer.

[Learn More](#)

Parameter	Value
Product Name	MGLL Rabbit pAb
Remarks/Alias	MGL; HUK5; MAGL; HU-K5; MGLL
Species	Human
Gene ID (Human)	11343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-313 of human MGLL (NP_009214.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99685
Protein Weight	33kDa
Shipping	Ice bag