


Anti-Mgmt Rabbit Monoclonal Antibody For Wb, Elisa - P16455

Item Number: CM0002305



Introduction

MGMT Rabbit mAb is a specific rabbit monoclonal antibody against human O-6-methylguanine-DNA methyltransferase, validated for Western blot and ELISA. Ideal for DNA repair and alkylation damage research.

[Learn More](#)

Parameter	Value
Product Name	MGMT Rabbit mAb
Remarks/Alias	MGMT; O-6-methylguanine-DNA methyltransferase
Species	Human
GeneID (Human)	4255
GeneID	4255
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 108-207 of human MGMT (NP_002403.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16455
Protein Weight	22kDa
Shipping	Ice bag