


Anti-Mgmt Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16455

Item Number: CM0007160



Introduction

MGMT Rabbit mAb for WB, IHC-P, IF/ICC, ELISA. Recognizes human O-6-methylguanine-DNA methyltransferase, a DNA repair enzyme. Suitable for cancer research and drug resistance studies.

[Learn More](#)

Parameter	Value
Product Name	MGMT Rabbit mAb
Remarks/Alias	MGMT; O-6-methylguanine-DNA methyltransferase
Species	Human
Gene ID	4255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-207 of human MGMT. (NP_002403.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P16455
Protein Weight	22kDa
Shipping	Ice bag