


Anti-Mgmt Rabbit Polyclonal Antibody For Wb, Elisa - P16455

Item Number: CM0017105



Introduction

Anti-MGMT rabbit polyclonal antibody suitable for WB and ELISA, detecting human and mouse MGMT. Ideal for DNA repair and drug resistance mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	MGMT Rabbit pAb
Remarks/Alias	MGMT
Species	Human
Gene ID (Human)	4255
Gene ID	4255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human MGMT (NP_002403.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16455
Protein Weight	22kDa
Shipping	Ice bag