

Anti-Mgp Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P08493

Item Number: CM0027027



Introduction

Rabbit polyclonal antibody targeting Matrix Gla protein (MGP), inhibitor of bone formation. Validated for WB, IHC-P, IF-P, ELISA; reacts with human, mouse, rat. Suitable for bone biology and calcification research.

[Learn More](#)

Parameter	Value
Product Name	MGP Rabbit pAb
Remarks/Alias	NTI; GIG36; MGLAP; MGP
Species	Human
GeneID (Human)	4256
GeneID	4256
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-103 of human MGP (NP_000891.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08493
Protein Weight	12kDa
Shipping	Ice bag