

Anti-Mgrn1 Polyclonal Antibody For Wb, If/Icc, Elisa - O60291

Item Number: CM0024028

Introduction

Rabbit polyclonal antibody targeting human MGRN1 (E3 ubiquitin-protein ligase) for WB, IF/ICC, ELISA. Validated cross-reactivity with mouse and rat. Ideal for endosome trafficking and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	MGRN1 Rabbit pAb
Remarks/Alias	RNF156; MGRN1
Species	Human
GeneID (Human)	23295
GeneID	23295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-538 of human MGRN1 (O60291).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60291
Protein Weight	61kDa
Shipping	Ice bag