

Anti-Mib1/Dip1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Yt6

Item Number: CM0025234



Introduction

Rabbit anti-MIB1/DIP1 polyclonal antibody validated for WB, ELISA; cross-reacts with mouse, rat. Targets E3 ubiquitin ligase in Notch pathway, apoptosis, ciliogenesis. 110 kDa.

[Learn More](#)

Parameter	Value
Product Name	MIB1/DIP1 Rabbit pAb
Remarks/Alias	MIB; DIP1; ZZZ6; DIP-1; LVNC7; ZZANK2; MIB1/DIP1
Species	Human
Gene ID (Human)	57534
Gene ID	57534
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 560-860 of human MIB1/DIP1 (NP_065825.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YT6
Protein Weight	110kDa
Shipping	Ice bag