

Anti-Mib2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Ax9

Item Number: CM0019274

Introduction

Rabbit polyclonal antibody targeting human MIB2, an E3 ubiquitin ligase in Notch signaling. Suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	MIB2 Rabbit pAb
Remarks/Alias	ZZZ5; ZZANK1; MIB2
Species	Human
GeneID (Human)	142678
GeneID	142678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-580 of human MIB2 (NP_001164158.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AX9
Protein Weight	104kDa
Shipping	Ice bag