

Anti-Mica Rabbit Polyclonal Antibody For Wb, Elisa - Q29983

Item Number: CM0013609



Introduction

Rabbit polyclonal antibody against human MICA (MIC-A, PERB11.1), validated for WB and ELISA. MICA is a stress-induced ligand for NKG2D, involved in immune surveillance. Shows cross-reactivity with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	MICA Rabbit pAb
Remarks/Alias	MIC-A; PERB11.1; MICA
Species	Human
GeneID (Human)	100507436
GeneID	100507436
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q29983
Protein Weight	43kDa
Shipping	Ice bag