


Anti-Micu1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bpx6

Item Number: CM0021743



Introduction

High-affinity Anti-MICU1 rabbit polyclonal antibody for WB and ELISA, reacting with human, mouse, and rat MICU1 (Q9BPX6). Ideal for mitochondrial calcium signaling and MCU complex research.

[Learn More](#)

Parameter	Value
Product Name	MICU1 Rabbit pAb
Remarks/Alias	CALC; EFHA3; MPXPS; CBARA1; ara CALC; MICU1
Species	Human
GeneID (human)	10367
GeneID	10367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-255 of human MICU1 (NP_001182447.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPX6
Protein Weight	54kDa
Shipping	Ice bag