

Anti-Micu2 Rabbit Polyclonal Antibody For Wb & Elisa - Q8Iyu8

Item Number: CM0013553



Introduction

High-quality rabbit polyclonal antibody targeting human MICU2 (EFHA1), mitochondrial calcium uptake regulator. Validated for WB and ELISA applications. Suitable for calcium signaling and mitochondrial function research.

[Learn More](#)

Parameter	Specification
Product Name	MICU2 Rabbit pAb
Remarks/Alias	EFHA1; 111008L20Rik; MICU2
Species	Human
GeneID (Human)	221154
GeneID	221154
Immunogen	Synthetic peptide. A synthetic peptide corresponding to a sequence within amino acids 200-300 of human MICU2 (NP_689939.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IYU8
Protein Weight	50kDa
Shipping	Ice bag