


Anti-Midkine (Mdk) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21741

Item Number: CM0018693



Introduction

Rabbit polyclonal antibody against human Midkine (MDK) for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for cytokine, growth factor, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	MDK Rabbit pAb
Remarks/Alias	MK; ARAP; NEGF2; MDK
Species	Human
GeneID (Human)	4192
GeneID	4192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-143 of human MDK (NP_001012333.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21741
Protein Weight	16kDa
Shipping	Ice bag