

Anti-Mief1 Rabbit Polyclonal Antibody (Mid51) For Wb And Elisa - Q9Nqg6

Item Number: CM0015003



Introduction

MIEF1 Rabbit pAb is a polyclonal antibody raised in rabbit against human MIEF1 (MID51), validated for Western blot (WB) and ELISA. It reacts with human and rat species. MIEF1 is a mitochondrial outer membrane protein that regulates fission/fusion, making this antibody ideal for mitochondrial dynamics research.

[Learn More](#)

Parameter	Value
Product Name	MIEF1 Rabbit pAb
Remarks/Alias	MID51, L0R8F8, SMCR7L, AltMIEF1, HSU79252, MIEF1-MP, dj1104E15.3, MIEF1
Species	Human
Gene ID (Human)	54471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-280 of human MIEF1 (NP_061881.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NQG6
Protein Weight	51kDa
Shipping	Ice bag