

Anti-Mien1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9BrT3

Item Number: CM0023206



Introduction

Rabbit polyclonal antibody targeting human MIEN1, validated for Western blot (WB) and ELISA. MIEN1 plays a role in cell migration and apoptosis, and is overexpressed in breast and prostate cancers. Ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	MIEN1 Rabbit pAb
Remarks/Alias	C35, ORB3, XTP4, RDX12, C17orf37, MIEN1
Species	Human
GeneID (Human)	84299
GeneID	84299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human MIEN1 (NP_115715.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BRT3
Protein Weight	12kDa
Shipping	Ice bag