

Anti-Mier1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N108

Item Number: CM0020109

Introduction

Rabbit polyclonal antibody against human MIER1, validated for Western blot, IHC-P, and ELISA. Recognizes human and mouse MIER1, a transcriptional repressor involved in chromatin silencing. Ideal for cancer research and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MIER1 Rabbit pAb
Remarks/Alias	ER1; MI-ER1; MIER1
Species	Human
GeneID (Human)	57708
GeneID	57708
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N108
Protein Weight	58kDa
Shipping	Ice bag