


Anti-Mif Monoclonal Antibody For Wb, Ihc-P, Elisa - P14174

Item Number: CM0007409



Introduction

High-quality rabbit monoclonal antibody against human MIF (Macrophage migration inhibitory factor). Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Plays a role in innate immune response and inflammation.

[Learn More](#)

Parameter	Value
Product Name	MIF Rabbit mAb
Remarks/Alias	GIF; GLIF; MMIF; MIF
Species	Human
GeneID (Human)	4282
GeneID	4282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-199 of human MIF (NP_002406.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14174
Protein Weight	12kDa
Shipping	Ice bag