

Anti-Miip Rabbit Polyclonal Antibody For Wb, Elisa - Q5Jxc2

Item Number: CM0018256

Introduction

Rabbit polyclonal antibody against human MIIP (IIP45) for Western blot and ELISA applications. Cross-reacts with rat. Suitable for cancer cell invasion and migration studies.

[Learn More](#)

Parameter	Specification
Product Name	MIIP Rabbit pAb
Remarks/Alias	IIP45; MIIP
Species	Human
Gene ID (Human)	60672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 279-388 of human MIIP (NP_068752.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q5JXC2
Protein Weight	43kDa
Shipping	Ice bag