

Anti-Mineralocorticoid Receptor (Nr3C2) Rabbit Pab - P08235

Item Number: CM0029293



Introduction

Rabbit polyclonal antibody against human Mineralocorticoid receptor (NR3C2), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. For mineralocorticoid signaling and hypertension research.

[Learn More](#)

Parameter	Value
Product Name	Mineralocorticoid receptor (NR3C2) Rabbit pAb
Remarks/Alias	MR; MCR; MLR; NR3C2VIT; Mineralocorticoid receptor (NR3C2)
Species	Human
Gene ID (Human)	4306
Gene ID	4306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-600 of human Mineralocorticoid receptor (NR3C2) (NP_000892.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08235
Protein Weight	107kDa
Shipping	Ice bag