


Anti-Mios (Mio, Sea4, Yulink) Polyclonal Antibody For Wb, Elisa - Q9Nxc5

Item Number: CM0023587



Introduction

Rabbit polyclonal antibody targeting MIOS (MIO, Sea4, Yulink), a GATOR2 complex component in mTOR signaling. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for cancer and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	MIOS Rabbit pAb
Remarks/Alias	MIO; Sea4; Yulink; MIOS
Species	Human
GeneID (Human)	54468
GeneID	54468
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 556-875 of human MIOS (NP_061878.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXC5
Protein Weight	99kDa
Shipping	Ice bag