


Anti-Miox Polyclonal Antibody For Wb, If-P, Elisa - Q9Ugb7

Item Number: CM0013233



Introduction

Rabbit polyclonal antibody against human MIOX (inositol oxygenase), validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for kidney-specific oxidoreductase studies.

[Learn More](#)

Parameter	Value
Product Name	MIOX Rabbit pAb
Remarks/Alias	ALDRL6; MIOX
Species	Human
GeneID (human)	55586
GeneID	55586
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human MIOX (NP_060054.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UGB7
Protein Weight	33kDa
Shipping	Ice bag