

Anti-Mip Polyclonal Antibody For Wb, Ihc-P, Elisa - P30301

Item Number: CM0029022

Introduction

Rabbit polyclonal antibody targeting human MIP (Aquaporin-0). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for lens water homeostasis and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	MIP Rabbit pAb
Remarks/Alias	AQP0; LIM1; MP26; MIP26; CTRCT15; MIP
Species	Human
GeneID (Human)	4284
GeneID	4284
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human MIP (NP_036196.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30301
Protein Weight	28kDa
Shipping	Ice bag