

Anti-Mipep Monoclonal Antibody For Wb, If/Icc, Elisa - Q99797

Item Number: CM0006547

Introduction

MIPEP rabbit monoclonal antibody validated for Western blot, IF/ICC, and ELISA. Targets mitochondrial intermediate peptidase (Q99797), ideal for studying mitochondrial protein processing and related diseases.

[Learn More](#)

Parameter	Value
Product Name	MIPEP Rabbit mAb
Remarks/Alias	MIP; HMIP; COXPD31; MIPEP
Species	Human
GeneID (Human)	4285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 511-713 of human MIPEP (NP_005923.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q99797
Protein Weight	81kDa
Shipping	Ice bag