

Anti-Mi pep Rabbit Polyclonal Antibody For Wb, Elisa - Q99797

Item Number: CM0013367

Introduction

MIPEP Rabbit pAb validated for WB and ELISA, cross-reactive with human, mouse, and rat. Targets mitochondrial intermediate peptidase (81 kDa), suitable for mitochondrial protein processing research.

[Learn More](#)

Parameter	Value
Product Name	MIPEP Rabbit pAb
Remarks/Alias	MIP; HMIP; COXPD31; MIPEP
Species	Human
GeneID (Human)	4285
GeneID	4285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 511-713 of human MIPEP (NP_005923.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99797
Protein Weight	81kDa
Shipping	Ice bag