

Anti-Mitd1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wv92

Item Number: CM0023949



Introduction

Rabbit polyclonal antibody against human MITD1, validated for WB, IHC-P, and ELISA. Detects a 29 kDa protein involved in cytokinesis and abscission. Suitable for human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	MITD1 Rabbit pAb
Remarks/Alias	MITD1
Species	Human
Gene ID (Human)	129531
Gene ID	129531
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WV92
Protein Weight	29kDa
Shipping	Ice bag