


Anti-Mitf Monoclonal Antibody For Wb, Ihc-P, Elisa - 075030

Item Number: CM0007243



Introduction

Rabbit monoclonal antibody targeting human MITF, validated for WB, IHC-P, and ELISA. MITF is a master transcription factor regulating melanocyte survival, differentiation, and pigmentation gene expression.

[Learn More](#)

Parameter	Value
Product Name	MITF Rabbit mAb
Remarks/Alias	MI; WS2; CMM8; WS2A; COMMAD; MITF-A; bHLHe32
Species	Human
GeneID (human)	4286
GeneID	4286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 402-526 of human MITF (NP_001341533.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	075030
Protein Weight	59kDa
Shipping	Ice bag