


Anti-Mitf Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O75030

Item Number: CM0007572



Introduction

Recombinant rabbit monoclonal antibody targeting human MITF (O75030), validated for WB, IF/ICC, ELISA. Reacts with human and rat MITF. Suitable for melanocyte differentiation, pigmentation gene regulation, and transcription factor research.

[Learn More](#)

Parameter	Value
Product Name	MITF Rabbit mAb
Remarks/Alias	MI; WS2; CMM8; WS2A; COMMAD; MITF-A; bHLHe32; MITF
Species	Human
Gene ID (Human)	4286
Gene ID	4286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 402-526 of human MITF (NP_001341533.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O75030
Protein Weight	59kDa
Shipping	Ice bag