


Anti-Mitf Rabbit Polyclonal Antibody For Wb And Elisa - 075030

Item Number: CM0013778



Introduction

Rabbit polyclonal antibody targeting human MITF (Microphthalmia-associated transcription factor). Validated for WB and ELISA, cross-reacts with human and mouse. MITF is a key regulator of melanocyte survival, differentiation, and melanosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	MITF Rabbit pAb
Remarks/Alias	MI; WS2; CMM8; WS2A; COMMAD; MITF-A; bHLHe32; MITF
Species	Human
Gene ID (Human)	4286
Gene ID	4286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human MITF (NP_937801.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75030
Protein Weight	59kDa
Shipping	Ice bag