

Anti-Mitofilin Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q16891

Item Number: CM0006513

Introduction

High-quality rabbit monoclonal antibody targeting human Mitofilin (MIC60) for Western blot, IHC, IP, and ELISA. Ideal for mitochondrial cristae morphology and MICOS complex research.

[Learn More](#)

Parameter	Value
Product Name	Mitofilin Rabbit mAb
Remarks/Alias	HMP; P87; P89; PIG4; Mic60; PIG52; MINOS2; P87/89; MICOS60; Mitofilin
Species	Human
GeneID	10989
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-300 of human Mitofilin (NP_006830.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16891
Protein Weight	84 kDa
Shipping	Ice bag