

Anti-Mitofusin-1/Mfn1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iwa4

Item Number: CM0025574



Introduction

Rabbit polyclonal antibody against human Mitofusin-1 (MFN1), a mitochondrial outer membrane GTPase essential for mitochondrial fusion. Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Suitable for mitochondrial dynamics research.

[Learn More](#)

Parameter	Value
Product Name	Mitofusin-1/MFN1 Rabbit pAb
Remarks/Alias	hfzo1; hfzo2; Mitofusin-1/MFN1
Species	Human
Gene ID (Human)	55669
Gene ID	55669
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 622-741 of human Mitofusin-1/MFN1 (NP_284941.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWA4
Protein Weight	84kDa
Shipping	Ice bag