

Anti-Mitofusin-1/Mfn1 Rabbit Monoclonal Antibody For Wb, IHC-P, Elisa - Q8Iwa4

Item Number: CM0004395



Introduction

High-quality rabbit monoclonal antibody against human Mitofusin-1 (MFN1), validated for WB, IHC-P, and ELISA. Detects 84 kDa protein, suitable for studying mitochondrial dynamics and fusion.

[Learn More](#)

Parameter	Value
Product Name	Mitofusin-1/MFN1 Rabbit mAb
Remarks/Alias	hfzo1; hfzo2; Mitofusin-1/MFN1
Species	Human
GeneID (human)	55669
GeneID	55669
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 642-741 of human Mitofusin-1/MFN1 (NP_284941.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q8IWA4
Protein Weight	84kDa
Shipping	Ice bag