

Anti-Mitofusin 2 (Mfn2) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95140

Item Number: CM0018424



Introduction

High-quality rabbit polyclonal antibody targeting human Mitofusin 2 (MFN2). Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for mitochondrial dynamics, fusion, and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Mitofusin 2 Rabbit pAb
Remarks/Alias	HSG; MARF; CMT2A; CPRP1; CMT2A2; HMSN6A; CMT2A2A; CMT2A2B; Mitofusin 2
Species	Human
Gene ID (Human)	9927
Gene ID	9927
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 70-260 of human Mitofusin 2 (NP_055689.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95140
Protein Weight	86kDa
Shipping	Ice bag