

Anti-Mki67 Rabbit Polyclonal Antibody For If-P, Ihc-P, Elisa - E9Pvx6

Item Number: CM0031077

Introduction

Rabbit polyclonal antibody against mouse Ki-67 (Mki67), validated for IF-P, IHC-P, and ELISA. Cross-reacts with rat. Ideal for cell proliferation and oncology research in IVD applications.

[Learn More](#)

Parameter	Value
Product Name	Mki67 Rabbit pAb
Remarks/Alias	Ki67; Ki-67; D630048A14Rik
Species	Mouse
GeneID (Human)	4288
GeneID	17345
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1816-2163 of mouse Mki67 (NP_001074586.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	E9PVX6
Protein Weight	350kDa
Shipping	Ice bag