

Anti-Mkks Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npj1

Item Number: CM0021403

Introduction

Rabbit polyclonal antibody against human MKKS (BBS6) for Western blot and ELISA. Recognizes a 62kDa molecular chaperone involved in BBSome assembly and ciliogenesis. Validated in human samples.

[Learn More](#)

Parameter	Specification
Product Name	MKKS Rabbit pAb
Remarks/Alias	KMS; MKS; BBS6; HMCS; MKKS
Species	Human
Gene ID (Human)	8195
Gene ID	8195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-570 of human MKKS (NP_740754.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NPJ1
Protein Weight	62kDa
Shipping	Ice bag