


Anti-Mkl1 Polyclonal Antibody For Wb, Elisa - Q969V6

Item Number: CM0018318



Introduction

Rabbit polyclonal antibody targeting human MKL1 (MRTF-A, MAL). Validated for WB and ELISA with cross-reactivity against mouse and rat. Ideal for studying cytoskeletal gene regulation, cell migration, and SRF-mediated transcription.

[Learn More](#)

Parameter	Value
Product Name	MKL1 Rabbit pAb
Remarks/Alias	MAL; MKL; BSAC; MKL1; MRTF-A
Species	Human
Gene ID (Human)	57591
Gene ID	57591
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human MKL1 (NP_065882.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969V6
Protein Weight	99kDa
Shipping	Ice bag