


Anti-Mkl1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q969V6

Item Number: CM0028678



Introduction

Rabbit polyclonal antibody targeting human MKL1 (MRTF-A/MAL/BSAC) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat MKL1. Ideal for cytoskeletal gene regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MKL1 Rabbit pAb
Remarks/Alias	MAL; MKL; BSAC; MKL1; MRTF-A
Species	Human
GeneID (Human)	57591
GeneID	57591
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human MKL1 (NP_065882.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969V6
Protein Weight	99kDa
Shipping	Ice bag