

Anti-Mkl2 (Mrtfb) Rabbit Polyclonal Antibody For Wb & Elisa - Q9ULh7

Item Number: CM0016867

Introduction

Rabbit polyclonal antibody raised against recombinant human MKL2 (MRTFB). Validated for Western blot (118 kDa) and ELISA, with cross-reactivity to human and mouse. Useful for research on SRF coactivation and muscle development.

[Learn More](#)

Parameter	Value
Product Name	MKL2 Rabbit pAb
Remarks/Alias	MKL2; MRTF-B; NPD001
Species	Human
GeneID (human)	57496
GeneID	57496
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 170-300 of human MKL2 (NP_054767.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULH7
Protein Weight	118kDa
Shipping	Ice bag