


Anti-Mkln1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UI63

Item Number: CM0018438



Introduction

Rabbit polyclonal antibody against human MKLN1 (Muskelin, TWA2), suited for WB and ELISA. Cross-reacts with mouse. Ideal for studying CTLH complex, cell spreading, and THBS1 signaling.

[Learn More](#)

Parameter	Value
Product Name	MKLN1 Rabbit pAb
Remarks/Alias	TWA2; MKLN1
Species	Human
Gene ID (Human)	4289
Gene ID	4289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 436-735 of human MKLN1 (NP_037387.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UL63
Protein Weight	85kDa
Shipping	Ice bag