


Anti-Mknk1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bub5

Item Number: CM0024320



Introduction

MKNK1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat MKNK1, a kinase regulating translation via EIF4E phosphorylation in stress and cytokine responses.

[Learn More](#)

Parameter	Value
Product Name	MKNK1 Rabbit pAb
Remarks/Alias	MNK1; MKNK1
Species	Human
Gene ID (Human)	8569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 374-465 of human MKNK1 (NP_003675.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUB5
Protein Weight	51 kDa
Shipping	Ice bag