

Anti-Mkp1 Rabbit Polyclonal Antibody For Wb, Elisa - P28562

Item Number: CM0023303



Introduction

Anti-MKP1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA. Detects human and mouse MKP1. Recombinant immunogen within amino acids 267-367. Ideal for studying MAP kinase signaling and stress response.

[Learn More](#)

Parameter	Value
Product Name	MKP1 Rabbit pAb
Remarks/Alias	HVH1; MKP1; CL100; MKP-1; PTPN10
Species	Human
Gene ID	1843
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 267-367 of human MKP1 (NP_004408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28562
Protein Weight	39kDa
Shipping	Ice bag