

Anti-Mkp1 Rabbit Polyclonal Antibody For Wb, Elisa - P28563

Item Number: CM0021745

Introduction

Rabbit polyclonal antibody targeting mouse MKP1 (DUSP1), validated for Western blot and ELISA. Reacts with human, mouse, and rat. MKP1 is a nuclear dual specificity phosphatase that dephosphorylates ERK2, regulating meiotic cell cycle.

[Learn More](#)

Parameter	Value
Product Name	MKP1 Rabbit pAb
Remarks/Alias	erp, MKP1, mkp-1, 3CH134, Ptpn16
Species	Mouse
GeneID	19252
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of mouse MKP1 (NP_038670.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28563
Protein Weight	39kDa
Shipping	Ice bag