

Anti-Mkrn3 Rabbit Polyclonal Antibody For Wb And Elisa - Q13064

Item Number: CM0023240

Introduction

MKRN3 Rabbit pAb is a polyclonal antibody derived from rabbit, recognizing human MKRN3 E3 ubiquitin ligase. Validated for Western blot and ELISA, it cross-reacts with human and mouse. For research on puberty timing and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	MKRN3 Rabbit pAb
Remarks/Alias	CPPB2; D1559; RNF63; ZFP127; ZNF127; MKRN3
Species	Human
GeneID (Human)	7681
GeneID	7681
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MKRN3 (NP_005655.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13064
Protein Weight	56kDa
Shipping	Ice bag