

Anti-Mks1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nxb0

Item Number: CM0030509

Introduction

MKS1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. 65 kDa protein. Ideal for cilia and ciliopathy studies. – Q9NXB0

[Learn More](#)

Parameter	Value
Product Name	MKS1 Rabbit pAb
Remarks/Alias	MES; MKS; BBS13; POC12; JBTS28; MKS1
Species	Human
GeneID (human)	54903
GeneID	54903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-300 of human MKS1 (NP_060247.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXB0
Protein Weight	65kDa
Shipping	Ice bag