


Anti-MLF1 Polyclonal Antibody For Wb, If/Icc, Elisa - P58340

Item Number: CM0028752



Introduction

Rabbit anti-MLF1 polyclonal antibody for WB, IF/ICC, ELISA, recognizing human and mouse MLF1. MLF1 regulates hematopoietic lineage commitment, cell cycle, and differentiation. Suitable for leukemia and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	MLF1 Rabbit pAb
Remarks/Alias	MLF1
Species	Human
GeneID (Human)	4291
GeneID	4291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human MLF1 (NP_071888.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P58340
Protein Weight	31kDa
Shipping	Ice bag