


## Anti-MLF1 Polyclonal Antibody For Wb, If/Icc, Elisa - P58340

Item Number: CM0013101



### Introduction

Rabbit polyclonal antibody against human MLF1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for research on myeloid leukemia factor 1 in hematopoietic differentiation, cell cycle regulation, and transcription factor function.

[Learn More](#)

| Parameter        | Specification                                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MLF1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | MLF1                                                                                                                                 |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 4291                                                                                                                                 |
| GeneID           | 4291                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human MLF1 (NP_071888.1). |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, IF/ICC, ELISA                                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                    |
| SWISS            | P58340                                                                                                                               |
| Protein Weight   | 31kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                              |