

Anti-MLF2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15773

Item Number: CM0016335

Introduction

Rabbit polyclonal antibody against human MLF2, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat MLF2. Suitable for leukemia and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	MLF2 Rabbit pAb
Remarks/Alias	NTN4; MLF2
Species	Human
GeneID (Human)	8079
GeneID	8079
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human MLF2 (NP_005430.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15773
Protein Weight	28kDa
Shipping	Ice bag