

Anti-MLh3 Monoclonal Antibody For Wb, Elisa - Q9Uhc1

Item Number: CM0004492

Introduction

Recombinant rabbit monoclonal antibody targeting human MLH3 (DNA mismatch repair protein Mlh3). Validated for WB and ELISA, with cross-reactivity to human and mouse. Ideal for DNA repair and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MLH3 Rabbit mAb
Remarks/Alias	HNPCC7
Species	Human
GeneID (Human)	27030
GeneID	27030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1190-1389 of human MLH3 (NP_001035197.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHC1
Protein Weight	164kDa
Shipping	Ice bag