

Anti-MLKL Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8Nb16

Item Number: CM0010885



Introduction

Mouse monoclonal antibody against human MLKL, validated for Western blot, IHC-P, and ELISA. Detects human and mouse MLKL, a 54 kDa necroptosis mediator. Ideal for cell death research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MLKL Mouse mAb
Remarks/Alias	hMLKL
Species	Human
GeneID (human)	197259
GeneID	197259
Immunogen	Synthetic peptide: A synthesized peptide derived from human MLKL.
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NB16
Protein Weight	54 kDa
Shipping	Ice bag