

Anti-MLKL Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q8Nb16

Item Number: CM0011895



Introduction

Rabbit monoclonal antibody against human MLKL (hMLKL) for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse. Ideal for necroptosis and programmed cell death research. Protein weight 54kDa.

[Learn More](#)

Parameter	Value
Product Name	MLKL Rabbit mAb
Remarks/Alias	hMLKL
Species	Human
GeneID (Human)	197259
GeneID	197259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 369-471 of human MLKL (NP_689862.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NB16
Protein Weight	54kDa
Shipping	Ice bag