


Anti-MLKL Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9D2Y4

Item Number: CM0026889



Introduction

Rabbit polyclonal antibody targeting mouse MLKL (mixed lineage kinase domain-like protein), validated for WB, IF-P, IHC-P, and ELISA applications. Cross-reacts with human MLKL. Ideal for necroptosis and cell death studies.

[Learn More](#)

Parameter	Value
Product Name	MLKL Rabbit pAb
Remarks/Alias	9130019I15Rik; MLKL
Species	Mouse
Gene ID (Human)	197259
Gene ID	74568
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-53 of mouse MLKL (Q9D2Y4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9D2Y4
Protein Weight	54kDa
Shipping	Ice bag